

Synaptosome and synaptosome fractionation

Rachel Jeffrey – adapted from Carlin RK, Grab DJ, Cohen RS, Siekevitz P (1980) Isolation and characterization of postsynaptic densities from various brain regions: enrichment of different types of postsynaptic densities. J Cell Biol 86:831-843 and MB Kennedy protocol (see any paper) and Mike Ehlers lab protocol (see Ehlerslab.org)

Preparations for the day before

Solution A (1X)	Final conc.	For 100 mL
Sucrose	0.32 M	10.944 g
NaHCO ₃ (0.5 M)	1 mM	200 uL
MgCl ₂ (1M)	1 mM	100 uL
CaCl ₂ ·2H ₂ O	0.5 mM	7.351 mg
Na pyrophosphate (FW: 265.9)	10 mM	265.9 mg
Nanopure water		
Protease inhibitors (add right before use)		2 tablets

Solution B (1X)	Final conc.	For 50 mL	For 100 mL
Sucrose	0.32 M	5.472 g	10.944 g
NaHCO ₃ (0.5 M)	1 mM	100 uL	200 uL
Nanopure water			
Protease inhibitors (add right before use)		1 tablet	

For gradient:

Concentration	1.2 M	1.0 M	0.85 M
Sucrose	8.208 g	6.844 g	5.836 g
1 mM NaHCO ₃	To 20 mL		

Note: sucrose adds a significant volume to the solution.

For PSD preps, also make 0.5% Triton X with Solution B

Note: Keep everything on ice, pre-cool reagents, add protease and phosphatase inhibitors to everything!

1. Rinse and dissect brains in **solution A** on ice
2. Homogenize with homogenizer using 12 up and down strokes. Use 20 mL **solution A** for 5 g brain tissue. One mouse forebrain is about 300 mg.
3. Combine homogenates and dilute to 10% weight/volume in **solution A** (about double of the original volume of solution A).
4. Centrifuge at low speed **710 g** in clinical centrifuge to get soft pellet for **10 minutes** (in 15 mL falcon tubes)
5. Resuspend pellet in the **solution A** supernatant using 3 up and down strokes of Teflon glass homogenizer.
6. Centrifuge at **1400 g** for **10 minutes** to remove nuclei.
7. Pool supernatants into 30 mL reusable centrifuge tubes and centrifuge **13,800 g** for **10 minutes** in floor centrifuge.
8. Resuspend pellet (containing synaptosomes and mitochondria) with six strokes of homogenizer in **solution B** (0.32 M sucrose, 1 mM NaHCO₃) using 16 mL solution B per 5 g starting material
9. Layer 8 mL of homogenized pellet on discontinuous **sucrose gradient** (10 mL in each layer of 1.2 M, 1.0 M, 0.85 M sucrose). Use 25 x 89 mm Beckman centrifuge tubes (reorder number 344058).
10. Spin in ultracentrifuge for **2 hours** at **82,500 g** (25,000 rpm for the SW28) at 4 °C. Set acceleration and brake to slow setting.
11. Aspirate top layers carefully with vacuum. Synaptosomes will be in layer between 1.0 M and 1.2 M; should get about ~3-5 mL.

To continue on to make PSDs:

1. Incubate synaptosomes in ice cold **0.5% Triton X in solution B** for **15 minutes** shaking on ice.
2. Centrifuge **35,000g** in Sorvall for **20 min** (PSD 1-Triton pellet).
3. Incubate this pellet again in 5 mL **Triton X solution** for **15 min** on ice.
4. Centrifuge **201,800g** for **1hr** to obtain PSD (two-Triton pellet) in SW41 rotor.
5. For protein extraction, solubilize pellet by drying completely, then add **SDS** up to final concentration of 3%. Sonicate pellet if necessary for **5-10 min**. Boil in **Sample buffer**.